NG-TEST® Blood Culture PREP

NG-TEST® Blood Culture PREP is an In Vitro Diagnostic (IVD) device reagent for sample preparation from blood culture, to be used in combination with NG-TEST® CARBA-5 and/or NG-TEST® CTX-M Multi.

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SK Návod na použitie - 52

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INTENDED USE

NG-TEST® Blood Culture PREP is an In Vitro Diagnostic (IVD) device reagent for sample preparation to be used in combination with **NG-TEST® CARBA-5** and/or **NG-TEST® CTX-M Multi**.

It is intended to prepare a bacterial sample from a positive blood culture with *Enterobacterales* or *Pseudomonas aeruginosa* from patients with a suspected BloodStream Infection (BSI).

The prepared sample is intended to be used in an IVD procedure to detect and differentiate carbapenemases and/or detect CTX-M groups extended-spectrum β -lactamases (ESBLs). The device is not automated and intended for professional use only.

SUMMARY

BSI is a life-threatening condition^{1,2} which is in most cases a propagation of an existing body-site infection to the blood. The aim of **NG-TEST® Blood Culture PREP** is to shorten the diagnosis procedure³ by preparing a bacterial sample to allow the detection of CTX-M ESBL-type or carbapenemases (KPC, OXA-48, VIM, IMP and NDM) directly from a positive blood culture without the need of subculture.

PRINCIPLE

NG-TEST® Blood Culture PREP contains a lysis buffer and a washing buffer. It is designed to prepare a bacterial sample ready to be used in combination with **NG-TEST® CARBA-5** and/or **NG-TEST® CTX-M Multi**.

The lysis buffer contains a detergent diluted in phosphate buffer with preservative and is mixed with the positive blood culture sample. This causes the release of hemoglobin (generating the red color of blood) from the red blood cells. The tube is centrifuged and the bacteria are recovered, free from hemoglobin and other blood components, into a pellet. The supernatant is discarded and the pellet is washed by adding washing buffer which contains a phosphate buffer with preservative. The tube is then centrifuged again and the supernatant discarded to recover a bacterial sample ready to be used with **NG-TEST® CTX-M Multi** or **NG-TEST® CARBA-5**.

REAGENTS AND MATERIALS SUPPLIED

NGP-BCP-S73-001

Each kit contains:

- 20 microtubes (2 mL).
- 1 Lysis buffer solution in a bottle with a red cap - R1 (23 mL).
- 1 Washing buffer solution in a bottle with a white cap - R2 (23 mL).
- 1 Instructions for use.

MATERIAL TO BE ORDERED SEPARATELY

NG-TEST® Blood Culture PREP is intended to be used in combination with **NG-TEST® CARBA-5** and/or **NG-TEST® CTX-M Multi**. The following references have to be order separately:

- NGB-CAR-S23-021 / **NG-TEST® CARBA-5**
- NGB-CAR-S23-022 / **NG-TEST® CARBA-5**
- NGB-CTM-S23-016 / **NG-TEST® CTX-M Multi**

MATERIAL REQUIRED BUT NOT SUPPLIED

- Centrifuge adapted to 2 mL microtubes
- Vortex
- Micropipette

PRECAUTIONS AND WARNINGS

- All reagents are for in vitro diagnostic use only.
- All the operations must be carried out in accordance with good laboratory practices.
- Wear protective gloves and safety glasses.
- Follow the procedure and precautions for use carefully to obtain accurate results. Do not mix up the steps during operating procedure.
- A blood-red final pellet is an indicator that the red cells have not been lysed/washed and indicates an error in the procedure.
- Do not use after the expiry date.
- Do not mix reagents from different kits.
- Close the buffer bottles tightly after use to prevent leakage. Make sure the cap is not screwed on at an angle.
- Inspect buffer bottles for damage or leaks before handling.

STORAGE AND STABILITY

Store the kit between 15°C and 30°C. The product is stable even after opening until the expiry date indicated on the label when stored in appropriate conditions.

Do not freeze and do not store outside the recommended temperature range.

Lysis buffer (solution R1 with red cap) may precipitate when stored at low temperatures. If precipitation occurs, see **Note 1** in operating procedure section.

SPECIMEN HANDLING AND COLLECTION

Samples must be obtained and handled in accordance with good microbiological practices.

Positive blood cultures must be collected in accordance with the manufacturer's recommendations.

NG-TEST® Blood Culture PREP can be used with the following blood culture bottles:

- bioMérieux: BACT/ALERT® FA Plus; BACT/ALERT® SA; BACT/ALERT® SN; BACT/ALERT® FN Plus; BACT/ALERT® PF Plus.
- Beckton Dickinson: BACTEC™ Plus Anaerobic medium; BACTEC™ Plus Aerobic medium; BACTEC™ Lytic Anaerobic medium; BACTEC™ Peds Plus medium

Positive blood cultures can be used up to 24h after being flagged positive either let under incubation or stored between 4°C and 30°C.

OPERATING PROCEDURE**PREPARATION :**

- Check the appearance of the Lysis buffer (solution R1 with red cap) to confirm the absence of precipitate.

NOTE 1: If precipitation occurs in Lysis Buffer R1, place the vial at 37°C for up to 1h with regular and gentle shaking until fully solubilized. Avoid strong shaking as this will create bubbles. Do not vortex. If the precipitate is still present, do not use the solution and contact your supplier.

- Under appropriate microbiological conditions, using the provided microtube, prepare an aliquot (>500 µL) of the positive blood culture.

REALISATION :

1. Gently homogenise the Lysis buffer (solution R1 with red cap) and add 1 mL of Lysis buffer into one of the provided microtube from the IVD test kit used in combination with. Avoid strong shaking as this will create bubbles. Do not vortex.
2. From the aliquot of positive blood culture, add 500 µL of the sample into the prepared microtube containing the Lysis buffer.
3. Close the microtube and vortex for at least 5 seconds.
4. Centrifuge for 1 minute, at 10,000 g (not refrigerated). Discard the supernatant by flipping the microtube upside down over a biological waste bin. Do not pipet the supernatant.

NOTE 2: After centrifugation, the bacterial pellet is attached to the bottom of the tube. It is not always visible.

5. Add 1 mL of washing buffer (solution R2 with white cap) into the microtube.
6. Close the microtube and vortex for at least 5 seconds to resuspend the pellet.
7. Centrifuge for 1 minute, at 10,000 g (not refrigerated). Discard the supernatant by flipping the microtube upside down over a biological waste bin. Do not pipet the supernatant.
8. Add the required number of drops of extraction buffer according to the IVD test used in combination and continue with the corresponding procedure.

PERFORMANCE

3 lots of **NG-TEST® Blood Culture PREP** were tested in combination with **NG-TEST® CARBA-5** and **NG-TEST® CTX-M Multi** to evaluate the inter-days, inter-lots and inter-operator's reproducibility and intra-runs repeatability with positive blood culture containing positive and negative strains for CTX-M enzymes or carbapenemases. Results are presented in the table below (percentage of compliant results).

TOTAL n=40	NG-TEST® CARBA-5		NG-TEST® CTX-M Multi	
SAMPLES	3 POSITIVES	1 NEGATIVE	2 POSITIVES	1 NEGATIVE
INTER-DAYS	100%	100%	100%	100%
INTER-OPERATORS	100%	100%	100%	100%
INTER-LOTS	100%	100%	100%	100%
INTRA-RUNS	100%	100%	100%	100%

LIMITATIONS

- The kit has only been validated for Gram negative **Enterobacterales** and **Pseudomonas aeruginosa**.
- Any use of the kit with another IVD procedure than the one indicated in the instructions for use should be validated first.
- The kit has only been validated with monomicrobial positive blood culture and should not be used with polymicrobial positive blood culture.

DISPOSAL

Dispose of all the materials used according to applicable local regulations regarding chemical and biohazardous waste.

Note 3: *Lysis buffer is classified as dangerous waste that must be disposed of in a special waste collection centre according the local regulations. Do not pour the buffer solutions in the sewer.*

TECHNICAL ASSISTANCE

Safety data sheets are available on request.

If you encounter a technical problem:

1. Record the lot number of the affected kit.
2. If possible, keep the sample and reagents in appropriate storage conditions while the complaint is processed.
3. Contact your local distributor or NG Biotech (support@ngbiotech.com).

Any serious incident that has occurred in relation to this device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

SAFETY INFORMATION

Lysis Buffer R1 is classified as a hazardous chemical: DANGER (H318) in accordance with EU Regulation N° 1272-2008 and EU Regulation N° 2015-830.

- H318 – Causes serious eye damage;
- P280 – Wear eye protection.
- P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 – Immediately call a POISON CENTER/doctor.

Contains sodium azid (<0.1 %)

BIBLIOGRAPHY

The bibliography can be found at the end of the document.

SYMBOLS

The symbols can be found at the end of the document.

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NG-BIOTECH® SAS