## NG•TEST®/CARBA-5

NG-TEST® CARBA-5 detects the presence of all the 3 classes of carbapenemase hydrolytic enzymes expressed by carbapenemase producing bacteria, identifying up to the 5 major types of carbapenemase enzymes in *Enterobacterales* or *Pseudomonas aeruginosa*, from a bacterial colony. The NG-TEST® CARBA-5 identifies the most prevalent carbapenemase families (NDM, IMP, VIM, OXA-48 and KPC) including their most clinically relevant variants, in the same cassette.

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## INTRODUCTION

**NG-TEST® CARBA-5** is a rapid in vitro diagnostic assay, for professional use only. It is a visual multiplex immunochromatographic (lateral flow) qualitative assay for the detection and differentiation of the five most common carbapenemase families (KPC, OXA-48-like, VIM, IMP and NDM) from carbapenemase producing, pure bacterial colonies of *Enterobacteriales* (including *Escherichia coli* and *Klebsiella pneumoniae*) and *Pseudomonas aeruginosa* in healthcare settings.

## SUMMARY AND PRINCIPLE

$\beta$ -Lactams are first-line antibiotics for the treatment of infections caused by *Enterobacteriales*; however, their efficacy has been challenged by the production of bacterial enzymes, known as  $\beta$ -lactamases, which hydrolyse the  $\beta$ -lactam ring and inactivate the antibiotic.  $\beta$ -Lactamases may be categorised into four classes according to Ambler's classification: classes A, C and D are  $\beta$ -lactamases with serine at the active site, whereas class B, also known as metallo- $\beta$ -lactamases or MBLs, have a zinc-dependent active site. Since the introduction of  $\beta$ -Lactam antibiotics into widespread clinical practice, the range of activity of  $\beta$ -Lactamases has grown.  $\beta$ -Lactamases with an extended-spectrum of activity (able to hydrolyse oxyimino-cephalosporins, and monobactams, but not cephamycins or carbapenems), known as extended-spectrum  $\beta$ -lactamases or ESBLs, emerged in hospital setting in the 1980s and became more common in community settings in the 1990s and early 2000s associated with the CTX-M family of ESBLs. Recently,  $\beta$ -lactamases with the ability to hydrolyze carbapenem antibiotics plus penicillins, cephalosporins and monobactam agents (known as carbapenemases) have been increasing in prevalence and have spread worldwide. Examples of the most prevalent carbapenemases found in *Enterobacteriales* are KPC (Class A), the metallo- $\beta$ -lactamases IMP, VIM and NDM (Class B), and OXA-48-like enzymes (Class D).

Together with clinical information and complementary diagnostic procedures, the use of **NG-TEST® CARBA-5** in the laboratory provides information that can support appropriate antimicrobial therapeutic choices, prompt infection control decision making and interventions.

**NG-TEST® CARBA-5** is an in vitro rapid and visual multiplex immunochromatographic assay that detects one or more of the five common types of carbapenemase enzymes (KPC (K), OXA-48-like (O), IMP (I), VIM (V), NDM (N)) from bacterial colonies. Extraction buffer is used as a cell lysing solution when mixed with colonies. Monoclonal antibodies that individually recognize each of the five carbapenemases are immobilized on a nitrocellulose membrane. Free monoclonal antibodies are present in the conjugate pad and labelled with colloidal gold. Upon addition of colonies mixed with extraction buffer to the sample pad, the capillary action of the nitrocellulose draws the sample through the mobile antibodies and immobile antibodies on the test strip. The immobilized control antibodies capture any mobile antibodies that run through the sample pad and nitrocellulose without binding to other test lines. A positive result occurs when a red line appears on the control region (C) and one or more lines appear in the test regions (K, O, V, I, or N) and indicates that the sample contains one or more carbapenemases. A negative result occurs when only the control line is observed and indicates that the sample does not contain any of the 5 carbapenemases. If the control line does not appear, the test result is invalid.

## REAGENTS AND MATERIALS SUPPLIED

### Each kit contains :

- 20 Test cassettes in aluminium pouches with desiccant
- 20 Microtubes
- 20 Disposable pipettes of 100  $\mu$ L
- 1 Extraction buffer solution in a plastic bottle (4,5 mL)

- 1 Instructions for use

### MATERIAL REQUIRED BUT NOT SUPPLIED

Timer • Single use gloves • Inoculation loop (1 µL) • Vortex

### PRECAUTIONS

- *In vitro* diagnostic test. For professional use only.
- All the operations must be carried out according to good laboratory practices.
- Do not use after the expiry date.
- The devices must remain in the sealed pouches until they are used.
- Handle the samples as if they were potentially infectious.
- After use, discard the device in an infectious waste container.
- Do not reuse the device.
- An inoculum size bigger than the one instructed in the operating mode is likely to cause false positive results on the test.

### STORAGE AND STABILITY

Store the devices in their sealed pouches between 4°C and 30°C. Do not freeze. Kits are stable in their intact packaging (pouch with desiccant) until the expiry date indicated on every kit. Buffer solution is reusable and stable when stored between 4°C and 30°C even after opening, until the expiry date indicated on the kit.

### CULTURE AND SAMPLING

The samples to be tested shall be obtained and handled according to standard microbiology procedures.

#### Validated culture media:

Luria Broth (LB) and LB agar, Trypticase soja agar (TSA), Mueller Hinton (MH) agar, UriSelect™ 4, Columbia agar + 5 % horse blood, Columbia agar + 5% sheep blood, ChromID® ESB� agar, ChromID® CARBA SMART, Drigalski (DRIG) agar, CHROMagar™ mSuperCARBA™, TSA + 5 % sheep blood, Mac Conkey, HardyCHROM™ CRE agar.

### OPERATING PROCEDURE

- 1 / Wear protective gloves.
- 2 / Bring the kit components at room temperature for at least 10 minutes.

#### Preparing the sample

- 1 / Dispense 5 drops of extraction buffer in one of the microtubes provided into the kit.
- 2 / From a solid agar-based culture, fill a 1 µL loop by touching on (not collecting) 3 colonies, and then suspend it in the microtube of extraction buffer.
- 3 / Close the microtube.
- 4 / Vortex to homogenise the mixture before use.

**NOTE:** *Viscous/mucous colonies can lead to migration problems, due to their high viscosity. When starting from viscous/mucous colonies, vortex for 3 minutes a colony in the extraction buffer and leave it to incubate for 10 minutes at room temperature before performing the test.*

#### Carrying out the test

- 1 / Open the pouch, and take out the device. Once opened, use the test immediately.
- 2 / Using the provided pipette, add 100 µL of the prepared mixture (sample must reach the black line indicated on the pipette to accurately dispense 100 µL) in the sample well labelled "S".
- 3 / Read the results at 15 minutes and interpret them as indicated below.

**NOTE:** *Do not interpret the test results after 20 minutes, as they may very possibly be causing false positive results.*

## RESULT INTERPRETATION


**Negative result**

If only one red line appears in the control region (C): the sample does not contain any of the indicated carbapenemase (or a non-detectable level) and must be interpreted as a negative result.

**Positive result**

If one red line appears in the control region (C) and one or several lines appear in the test regions K, O, V, I, N: the sample contains one or several carbapenemases and must be interpreted as a positive result.

**NOTE:** The intensity of the red test line(s) may vary depending on the concentration of the carbapenemases present in the sample. A weak line should be considered as a positive result.

**Invalid result**

If the control line (C) does not appear, or it is not aligned with the letter (C), the test result is invalid. Insufficient sample volume or an incorrect procedure are the most likely reasons for control line failure.

Deterioration of the test kit may have occurred. Repeat the procedure using a new test. If the problem persists, do not reuse the kit and contact your distributor.

## QUALITY CONTROL

An internal quality control is included in the test. When the control line develops, it confirms the sample volume was sufficient and the procedure was correct.

## LIMITATIONS

- 1 / A negative result does not preclude the presence of carbapenemase producing organisms (example: producing SME, GES, IMI).
- 2 / False negative results may occur with multiple subcultures of a bacterial isolate without any selective pressure.
- 3 / This test has not been validated for *Acinetobacter*.
- 4 / Some false negative results may occur with *Proteus mirabilis* tested from blood agar only.
- 5 / *Proteus* spp. tend to exhibit swarming growth on blood agar, thus only the surface of the swarming should be touched with a loop in 3 different colonies.
- 6 / This test is a qualitative assay and will not yield any quantitative result.
- 7 / This test should be used as an aid for the rapid identification of patients bearing a resistance to carbapenem antibiotics. The obtained results must be confirmed with alternative or complementary diagnostic procedures.
- 8 / A positive or a negative test does not rule out the presence of other mechanisms of antibiotic resistance.

## PERFORMANCES AND CHARACTERISTICS

**Detection limit:**

The detection limit was determined using purified recombinant enzymes:  
KPC 600 pg/mL • OXA 300 pg/mL • VIM 300 pg/mL • IMP 200 pg/mL • NDM 150 pg/mL

**Clinical Evaluation:**

**NG-TEST® CARBA-5** was evaluated at the NRC (National Reference Center for Antibiotic Resistance, Kremlin-Bicêtre Hospital, Paris, France). 147 strains were blind-tested with WGS-characterized  $\beta$ -lactamase content. Among the evaluated collection, only 2 rare D179Y KPC variants were not detected (KPC-31 and KPC-33). This D179Y impaired the carbapenem hydrolytic activity of the enzyme which cannot be considered as a carbapenemase, leading also to negative result by biochemical tests. The sensitivity and specificity of the **NG-TEST® CARBA-5** was calculated by considering the KPC-31 and KPC-33 variants as negative results.

**Table 1: Results obtained at NRC Kremlin Bicêtre<sup>11</sup>**

		WGS Whole Genome Sequencing		
		POSITIVE	NEGATIVE	TOTAL
NG•TEST® CARBA-5	POSITIVE	119	0	119
	NEGATIVE	0	28	28
	TOTAL	119	28	147
SENSITIVITY		100%	CI 95% = 96.9 - 100%	
SPECIFICITY		100%	CI 95% = 87.9 - 100%	

**Summary of variants detected by NG-TEST® CARBA-5 during validations:**

Organism group	Target	Variants detected during validations 8; 11; 12
<i>Enterobacteriales*</i>	KPC	2,3,4,5,6,7,12,14,28,39
	OXA-48-like	48, 162, 163, 181, 204, 232, 244, 245, 370, 405, 436, 484, 505, 515, 517, 519, 535, 793
	VIM	1, 2, 4, 5, 6, 11, 19, 23, 27, 31
	IMP	1, 4, 6, 8, 10, 11, 14, 26, 47, 58
	NDM	1, 3, 4, 5, 6, 7, 9, 11, 19
<i>Pseudomonas aeruginosa</i>	KPC	2, 5
	OXA-48-like	181
	VIM	2, 5, 11
	IMP	1, 2, 5, 7, 13, 14, 15, 16, 18, 19, 26, 29, 31, 37, 39, 46, 56, 63, 71, 79
	NDM	1

## Summary of variants detected by NG-TEST® CARBA 5 in publications:

Organism group	Target	Variants detected in publications 8; 9; 10; 11; 12; 13; 14; 15
<i>Enterobacteriales*</i>	KPC	2, 3, 4, 5, 6, 7, 12, 14, 19, 23, 25, 28, 35, 39, 51, 78, 79, 82, 93
	OXA-48-like	48, 162, 163, 181, 204, 232, 244, 245, 370, 405, 436, 484, 505, 515, 517, 519, 535, 793
	VIM	1, 2, 4, 5, 6, 11, 19, 23, 26, 27, 31, 39, 46, 51, 52, 54, 56, 58, 59, 67
	IMP	1, 4, 6, 7, 8, 10, 11, 14, 22, 26, 30, 47, 58, 61
	NDM	1, 2, 3, 4, 5, 6, 7, 8, 9, 11, 19
<i>Pseudomonas aeruginosa</i>	KPC	1, 2, 5
	OXA-48-like	181
	VIM	2, 5, 11
	IMP	1, 2, 4, 5, 7, 8, 13, 14, 15, 16, 18, 19, 25, 26, 29, 31, 37, 39, 45, 46, 56, 62, 63, 71, 79
	NDM	1

\*Including some: *Citrobacter braakii*, *C. freundii*, *C. koseri*, *C. sedlakii*, *Enterobacter aerogenes*, *E. asburiae*, *E. cloacae*, *E. cloacae* complex, ***Escherichia coli***, *E. hermannii*, *Klebsiella aerogenes*, *K. variicola*, *K. ozaenae*, *K. oxytoca*, ***K. pneumoniae***, *Kluyvera ascorbate*, *Morganella morganii*, *Proteus mirabilis*, *Providencia stuartii*, *P. rettgeri*, *Pantoea* species, *Providencia alcalifaciens*, *Raoultella ornithinolytica*, *Salmonella enterica*, *S. senftenberg*, *Serratia marcescens*, *Shigella boydii*, *S. sonnei*, *S. bicestris*, *Yersinia enterocolitica*.

Non-detection of KPC-31, KPC-33, KPC-52, KPC-71, KPC-76, KPC-77 variants.

## BIBLIOGRAPHY

The bibliography can be found at the end of the document.

## SYMBOLS

Content for  
20 assays

*in vitro* diagnostic  
medical device

Batch number

Importer

Expiry date

Do not re-use

Catalogue  
reference

Manufacturer

Temperature  
limit

Consult  
instructions  
for use

Unique Device  
Identifier

This kit was developed in collaboration with the  
CEA\*.

\*The French Alternative Energies and Atomic  
Energy Commission is a key player in research,  
development and innovation.

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1. Yohei D et al. Carbapenemase-Producing Enterobacteriaceae. *Semin Respir Crit Care Med.* 2015; 36(1): 74–84.
2. Nordmann P et al. Rapid Detection of Carbapenemase-producing Enterobacteriaceae. *Emerging Infectious Diseases.* 2012; 18(9).
3. Bush K et al. Updated Functional Classification of  $\beta$ -Lactamases. *Antimicrobial Agents and Chemotherapy.* 2010; 54(3):969-976.
4. Nordmann P et al. Carbapenemases in enterobacteriaceae. *Archives de Pédiatrie.* 2010;17:S154-S162.
5. Demir Y et al. Investigation of VIM, IMP, NDM-1, KPC AND OXA-48 enzymes in Enterobacteriaceae strains. *Pak J Pharm Sci.* 2015; May;28(3 Suppl):1127-33.
6. Oteo J et al. Evolution of carbapenemase-producing Enterobacteriaceae at the global and national level: what should be expected in the future? *Enferm Infecc Microbiol Clin.* 2014; 32 Suppl 4:17-23.
7. Kulkarni MV et al. Use of imipenem to detect KPC, NDM, OXA, IMP, and VIM carbapenemase activity from gram-negative rods in 75 minutes using liquid chromatography-tandem mass spectrometry. *J Clin Microbiol.* 2014; 52(7):2500-5.
8. Boutal H et al. 2018. A multiplex lateral flow immunoassay for the rapid identification of NDM-, KPC-, IMP- and VIM-type and OXA-48-like carbapenemase-producing Enterobacteriaceae. *J Antimicrob Chemother.*
9. Baeza LL, et al. A Comparison of five methods for detection of carbapenemases in Entérobactéries with proposal of a new algorithm, *Clinical Microbiology and Infection,*
10. Hopkins KL, Meunier D, Naas T, Volland H, Woodford N. 2018. Evaluation of the NG-TEST® CARBA-5 multiplex immunochromatographic assay for the detection of KPC, OXA-48-like, NDM, VIM and IMP carbapenemases. *J Antimicrob Chemother* 73:3523- 3526.
11. Volland H. et al., Improvement of the immunochromatographic NG-TEST® CARBA-5 assay for the detection of IMP-variants previously undetected *Antimicrobial Agents and Chemotherapy.* 2019
12. Klavins A et al. Evaluation of NG-TEST® CARBA-5 with Carbapenemase producing organisms and non-carbapenemase producing organisms. CPHM 991. *ASM Microbe* 2019-Chicago.
13. Hazim O. et al. Comparative Evaluation of Five Assays for Detection of Carbapenemases with a Proposed Scheme for Their Precise Application, *The Journal of Molecular Diagnostics* 2020; 22 (9)
14. Gu D, Yan Z, Cai C, Li J, Zhang Y, Wu Y, Yang J, Huang Y, Zhang R, Wu Y. Comparison of the NG-Test Carba 5, Colloidal Gold Immunoassay (CGI) Test, and Xpert Carba-R for the Rapid Detection of Carbapenemases in Carbapenemase-Producing Organisms. *Antibiotics (Basel).* 2023 Feb 2;12(2):300
15. Mendez-Sotelo BJ, López-Jácome LE, Colín-Castro CA, Hernández-Durán M, Martínez-Zavaleta MG, Rivera-Buendía F, Velázquez-Acosta C, Rodríguez-Zulueta AP, Morfín-Otero MDR, Franco-Cendejas R. Comparison of Lateral Flow Immunochromatography and Phenotypic Assays to PCR for the Detection of Carbapenemase-Producing Gram-Negative Bacteria, a Multicenter Experience in Mexico. *Antibiotics (Basel).* 2023 Jan 6;12(1):96.

# Importer information

## Australia

Distributor: Cell biosciences  
Pty Ltd  
PO Box 243, Heidelberg, VIC 3084  
Fax: 03 8678 3956  
Mobile: 0394 167 177  
Web: [www.cellbiosciences.com.au](http://www.cellbiosciences.com.au)

Note: any serious incident that occurs in relation to the device should be reported to the manufacturer and to the Therapeutic Goods Agency (TGA).

TGA:  
TGA, PO Box 100, Woden ACT 2606,  
Australia.  
Website: <https://www.tga.gov.au/>  
Email: [info@tga.gov.au](mailto:info@tga.gov.au)  
Fax: 02 6203 1605

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**NG·BIOTECH® SAS**