

NG-Test[®] Reader[®]



Patent pending



Features :

- ☛ Mobility
- ☛ Connectivity
- ☛ Traceability
- ☛ Data management

Applications :

- ☛ Human Health
- ☛ Animal Health
- ☛ Food Industry
- ☛ Bio Defense

NG-Test[®] Reader[®] Reliable tool for reading and providing quantitative results for your lateral flow rapid tests.



General procedure

- 1 Rapid test is inserted into Reader
- 2 Reader communicates image to a Mobile App (Bluetooth) or PC (USB)
- 3 Image is analyzed and the test signals are translated into a result
- 4 Result is managed by Mobile App or PC Software User Interface
- 5 Results can be communicated enabling digital data management